



Email : sales@laboratoryinstrumentindia.com

Product Name :
PIP Uniform Technology Reusable (Class 10) Cleanroom
Coverall

Product Code :
LBNY-00017-11000607



Description :

Coverall that Groundable Static Control Garment.
Coil zipper
Snaps at top of zipper
Raglan sleeves for more movement
Knit cuff
Color coded hang loop
Snaps at leg hem.

Technical Specification :

PIP Uniform Technology Reusable (Class 10) Cleanroom Coverall



Laboratory Instrument India