



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
PIP Uniform Technology Reusable (Class 10) Cleanroom  
Coverall

**Product Code :**  
LBNY-00017-11000607



**Description :**

Coverall that Groundable Static Control Garment.  
Coil zipper  
Snaps at top of zipper  
Raglan sleeves for more movement  
Knit cuff  
Color coded hang loop  
Snaps at leg hem.

**Technical Specification :**

PIP Uniform Technology Reusable (Class 10) Cleanroom Coverall



**Laboratory Instrument India**