



Email : sales@laboratoryinstrumentindia.com

Product Name :
Protective Products ATA Cut Resistant Sleeve 2-ply Design

Product Code :
LBNY-00017-11000616



Description :

Sta-Cool inner layer wicks moisture away from the arm, keeping the wearer dry, cool and comfortable. Keep the wrist protected.

Durability, strong cut and abrasion resistance are provided by ATA fiber technology

Sta-Cool inner layer wicks moisture away from the arm, keeping the wearer dry, cool and comfortable all day

Serious Arm Protection, All Day Comfort

Keep arms cool, comfortable and well protected with this ANSI A4 cut resistant sleeve

Machine washable.

Outstanding stretch means this sleeve will conform to just about any size arm – big or small – and can be worn comfortably over clothing.

Technical Specification :

Protective Products ATA Cut Resistant Sleeve 2-ply Design



Laboratory Instrument India