



Email : sales@laboratoryinstrumentindia.com

Product Name :
West Chester White Nylon Hairnets

Product Code :
LBNY-00017-11000618



Description :

Ideal for use in food service, food prep, manufacturing, and retail applications
Elastic band for flexible fit
Nylon hairnets feature an elastic adjustment and is available in 18 in. or 21 in. diameter.
Nylon hairnet
Available in 18 in. or 21 in. sizes
White.

Technical Specification :

West Chester White Nylon Hairnets



Laboratory Instrument India