



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Medium Weight Terrycloth Sleeves

**Product Code :**  
LBNY-00017-11000619



**Description :**

Medium Weight Terrycloth Sleeves offer heat resistant protection for up to 200°F.  
Elastic top/bottom provide secure fit and protects against harmful particles from entering the sleeve.  
Medium weight construction provides great comfort and flexibility while providing heat protection to the arm.  
Elastic inserted at top and bottom.  
Protects up to 200°F  
White seamless sleeve, loop-in.

**Technical Specification :**

Medium Weight Terrycloth Sleeves



**Laboratory Instrument India**