



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Plug Key

**Product Code :**  
LBNY-0003-1010021



**Description :**

Plug Key

**Technical Specification :**

Along with accurately tapered brass plugs with a fluted plastic top to fit in to the hole between the blocks, thus providing negligible contact resistance.  
Comprises two brass blocks 8mm thick fitted with brass 4mm socket terminals mounted in a moulded base with dimension 100 x 63 x 43mm



**Laboratory Instrument India**