



Email : sales@laboratoryinstrumentindia.com

Product Name :
High-Tech Conversions Level 2 Isolation Gown, Blue

Product Code :
LBNY-00017-11000624



Description :

Level 2 Isolation Gown is ideal for a variety of health care settings where fluid and splash protection is needed. Made from laminated spun bond polypropylene, these gowns are fluid resistant, full length, knit cuff, latex free and disposable.

Technical Specification :

High-Tech Conversions Level 2 Isolation Gown, Blue



Laboratory Instrument India