



Email : sales@laboratoryinstrumentindia.com

Product Name :
High-Tech Conversions Bee-Safe Open Face Cleanroom Hood

Product Code :
LBNY-00017-11000625



Description :

Provides complete head coverage and is made of lightweight, low-linting, and breathable polypropylene material. Open face cleanroom hood is lightweight with particulate filtration. Polypropylene cleanroom hood is ideal for a variety of applications including light industrial, lab and cleanroom. Latex-free.

Technical Specification :

High-Tech Conversions Bee-Safe Open Face Cleanroom Hood



Laboratory Instrument India