



Email : sales@laboratoryinstrumentindia.com

Product Name :
High-Tech Conversions Bee-Safe Cleanroom
Polypropylene Boot Covers

Product Code :
LBNY-00017-11000626



Description :

Polypropylene boot covers feature PE outer layer to protect against outside contaminants.
Polyethylene outer layer will ensure protection from outside contaminants
Lightweight material allows for longer wear
Manufactured from 100% polypropylene material providing a soft, comfortable and low-linting solution
Latex-free elastic band allows for a secure fit and closure
Super sticky PVC sole reduces slipping and sliding
Production process guarantees extremely low particle generation and chemical extractables.

Technical Specification :

High-Tech Conversions Bee-Safe Cleanroom Polypropylene Boot Covers



Laboratory Instrument India