



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
VAI 1700 MP Cleanroom Lab Coats

**Product Code :**  
LBNY-00017-11000630



**Description :**

MP lab coats are constructed of a non-woven, microporous, coated material and provide protection in controlled environments.

Liquid chemical resistance

Comfortable, yet durable

Bonded polypropylene with laminate

Low particle shedding

Disposable.

**Technical Specification :**

VAI 1700 MP Cleanroom Lab Coats



**Laboratory Instrument India**