



Email : sales@laboratoryinstrumentindia.com

Product Name :

VAI Easy2Gown 1700 Cleanroom Hoods

Product Code :

LBNY-00017-11000631



Description :

Folded in an inside-out fashion that essentially eliminates the individual from contaminating the outside of the garment

VAI Easy2Gown 1700 Cleanroom Hoods are single-use PPE with excellent chemical resistance

Bonded polypropylene with laminate

Piped seams with reinforced stitching

Packaged to prevent operator contamination of the outside of the garment

Excellent chemical resistance.

Technical Specification :

VAI Easy2Gown 1700 Cleanroom Hoods



Laboratory Instrument India