



Email : sales@laboratoryinstrumentindia.com

Product Name :

VAI Sterile MP Hood, Integrated Face Mask

Product Code :

LBNY-00017-11000632



Description :

Heavy weight, non-woven and high quality microporous (MP) coated material
Suited for use where maximum protection to product or wearer is required
Microporous hood with integrated face mask is individually packaged and suitable for ISO 5 cleanroom environments.
Liquid resistance, especially to chemicals.

Technical Specification :

VAI Sterile MP Hood, Integrated Face Mask



Laboratory Instrument India