



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
High-Tech Conversions Bee-Safe Level 1 Isolation Gowns

**Product Code :**  
LBNY-00017-11000642



**Description :**

Made from Spunbonded Polypropylene, these gowns are fluid resistant, full length, knit cuff, latex free and disposable.

Bee-Safe Level 1 Isolation Gowns are ideal for a variety of health care settings.

Material: Spunbonded Polypropylene

Sizes Available: Large (GISO1YL), X-Large (GISO1YX)

Additional Details: No pockets.

**Technical Specification :**

High-Tech Conversions Bee-Safe Level 1 Isolation Gowns



**Laboratory Instrument India**