



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
MP Protective Sleeves

**Product Code :**  
LBNY-00017-11000644



**Description :**

Heavyweight, non-woven and high quality microporous (MP) coated material  
Heavyweight, non-woven, microporous material is low in particulate shedding and suited for use where maximum protection is needed.  
Suited for use where maximum protection to product or wearer is required  
Liquid and chemical resistant.

**Technical Specification :**

MP Protective Sleeves



**Laboratory Instrument India**