



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

VAI Easy2Gown 1700 SMS Classic Cleanroom Coveralls

**Product Code :**

LBNY-00017-11000647



**Description :**

Breathable 1700 coveralls are made of a heavyweight, non-woven, and high-quality Micro Porous coated material. Disposable garments are extremely tough and fluid resistant.

Piped seams with reinforced stitches

Excellent chemical resistance

Bonded polypropylene with laminate

Sewn thumb loops to eliminate shifting

Disposable.

**Technical Specification :**

VAI Easy2Gown 1700 SMS Classic Cleanroom Coveralls



**Laboratory Instrument India**