



Email : sales@laboratoryinstrumentindia.com

Product Name :
Brady TIME spot Badge

Product Code :
LBNY-00017-11000649



Description :

Brady TIME spot badges have a self-expiring circle that lasts one day, preventing repeated use of a badge. Choose Visitor, Contractor or Temporary badges. Minimize the unauthorized use of badges. Apply the Time token or TIME spot badge over the printed "back" to start expiration. Low-cost clip-on paper badges with time-dependent "backs" printed on them. Badges expire in one day. 500 per package.

Technical Specification :

Brady TIME spot Badge



Laboratory Instrument India