



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Single Fixed Pulley

**Product Code :**  
LBNY-0003-1010024



**Description :**

Single Fixed Pulley

**Technical Specification :**

Made of Aluminium.  
Load voltage: \$ number kn static, 10kn dynamic.  
Fit for maximum rope diameter: 13mm.  
Half circle point of the hand, forced equally.  
Wheel pulley spillage, smooth the polishing surface.



**Laboratory Instrument India**