



Email : sales@laboratoryinstrumentindia.com

Product Name :
Premium Brilliant Ladies Fitted Vest

Product Code :
LBNY-00017-11000666



Description :

Features extra reflective piping, form fittedness for ladies.
High performance reflective material (Width: 2 in.)
Trim, reflective piping around arms, down sides for added low light visibility
Ultra-Cool 100% polyester mesh construction
Zipper front closure
Waist cinches to ensure perfect fit
Back, sides with lightweight solid front.

Pockets:

One inside ergonomically positioned lower left radio pocket with flap
Two hidden lower outside slash pockets with zipper closures.

Technical Specification :

Premium Brilliant Ladies Fitted Vest



Laboratory Instrument India

