



Email : sales@laboratoryinstrumentindia.com

Product Name :
Premium Brilliant Series Heavy Duty Vest

Product Code :
LBNY-00017-11000668



Description :

Made of Ultra-Cool 100% polyester mesh
Brilliant Trim reflective edging improves low light visibility
Zipper front closure reinforced with durable webbing
Created for maximum visibility. ML Premium Brilliant Series Heavy Duty Vest enhances low light visibility and provides instant differentiation between worker and the work zone.
Features left and right chest mic tabs
2 in. wide high performance reflective material with 3 in. breathable contrasting color

Pockets:

Two lower outside cargo pockets with grommets and adjustable flaps
Two inside patch pockets.
Right chest 2-tier pencil pocket
Left chest radio pocket with flap.

Technical Specification :

Premium Brilliant Series Heavy Duty Vest

Laboratory Instrument India

