



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Premium Brilliant Series Heavy Duty Vest

**Product Code :**  
LBNY-00017-11000668



**Description :**

Made of Ultra-Cool 100% polyester mesh  
Brilliant Trim reflective edging improves low light visibility  
Zipper front closure reinforced with durable webbing  
Created for maximum visibility. ML Premium Brilliant Series Heavy Duty Vest enhances low light visibility and provides instant differentiation between worker and the work zone.  
Features left and right chest mic tabs  
2 in. wide high performance reflective material with 3 in. breathable contrasting color

**Pockets:**

Two lower outside cargo pockets with grommets and adjustable flaps  
Two inside patch pockets.  
Right chest 2-tier pencil pocket  
Left chest radio pocket with flap.

**Technical Specification :**

Premium Brilliant Series Heavy Duty Vest

---

# Laboratory Instrument India

