



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Slinky Spring

**Product Code :**  
LBNY-0003-1010026



**Description :**

Slinky Spring

**Technical Specification :**

To illustrate wave motion. Flat section steel wire in helical coil.

The Coil Slinky spring is perfect for classroom demonstrations of wave properties. Great for many experiments, hands-on, a student favorite

Perfect for stocking stuffers, prizes, presents, and more because everyone loves a Slinky.

Coiled, diameter: 78mm length, closed 115mm.

Fully extended 5m.



**Laboratory Instrument India**