



Email : sales@laboratoryinstrumentindia.com

Product Name :
Slinky Spring

Product Code :
LBNY-0003-1010026



Description :

Slinky Spring

Technical Specification :

To illustrate wave motion. Flat section steel wire in helical coil.
The Coil Slinky spring is perfect for classroom demonstrations of wave properties. Great for many experiments, hands-on, a student favorite
Perfect for stocking stuffers, prizes, presents, and more because everyone loves a Slinky.
Coiled, diameter: 78mm length, closed 115mm.
Fully extended 5m.



Laboratory Instrument India