



Email : sales@laboratoryinstrumentindia.com

Product Name :
FR Lab Coat with Knit Cuffs

Product Code :
LBNY-00017-11000676



Description :

This lab coat is designed for functionality, comfort and safety with a CAT 1 arc rating and is also compliant.
Four front mitered style pockets for ample storage
Pass-through at hips to access underneath
Wash Care: Industrial Laundry - Light Soil
Rib cuffs for protection during movement
Concealed snap front closure to avoid scratching surfaces
Hazard: Arc Flash, Flash Fire
Protection: ATPV Arc Rating: 4.4.

Technical Specification :

FR Lab Coat with Knit Cuffs



Laboratory Instrument India