



Email : sales@laboratoryinstrumentindia.com

Product Name :
Solar Training Set (Power Case)

Product Code :
LBNY-0002-620026



Description :

Solar Training Set (Power Case)

Technical Specification :

Solar Energy Trainer is a versatile training system to be used in laboratories. It introduces the basic concept of solar cell (photovoltaic cell) of converting sunlight directly into electrical energy. The system works as a live demonstration unit of solar cells in terms of its concepts and applications.

This system is provided with Solar Energy Trainer and a Solar Panel.

Solar Training Set
Solar Plate & Inverter Drive

Solar Training Set (Power Case), Solar Training Set (Power Case) Equipments, Solar Training Set (Power Case) Tools, Solar Training Set (Power Case) Tool Kits, Solar Training Set (Power Case) Manufacturers, Solar Training Set (Power Case) Suppliers from India, China, Kenya



Laboratory Instrument India

