



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Battery Holder

**Product Code :**  
LBNY-0003-1010027



**Description :**

Battery Holder

**Technical Specification :**

PC Mountable and wire lead

High Performance Leaf Spring Holders

It can hold standard size batteries, either for alkaline or rechargeable batteries, the pin design is great for easily soldering and connecting.

Holder Case DIY Lithium Battery Box Battery Holder with Pin PCB Mounting type.

Plastic box, 10 of two cells and 10 of one cell.



**Laboratory Instrument India**