



Email : sales@laboratoryinstrumentindia.com

Product Name :
Battery Holder

Product Code :
LBNY-0003-1010027



Description :

Battery Holder

Technical Specification :

PC Mountable and wire lead
High Performance Leaf Spring Holders
It can hold standard size batteries, either for alkaline or rechargeable batteries, the pin design is great for easily soldering and connecting.
Holder Case DIY Lithium Battery Box Battery Holder with Pin PCB Mounting type.
Plastic box, 10 of two cells and 10 of one cell.



Laboratory Instrument India