



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
International Enviroguard Micro Guard MP Microporous  
Shoe Cover with Elastic Closure

**Product Code :**  
LBNY-00017-11000683



**Description :**

Lightweight, tear-resistant liquid and particulate barrier helps safeguard from blood penetration, various industrial chemicals, along with dust, dirt, and grime.

Using state-of-the-art lamination technology, Micro Guard MP microporous shoe covers combine a soft, non-woven fabric with a microporous film.

Shoe covers come with a strong elastic closure at the top to keep them in place.

Breathable fabric allows moisture to escape, so your feet are more cool and comfortable even as you sweat.

Environmentally friendly; can be recycled after use.

Passes testing for blood penetration.

**Technical Specification :**

International Enviroguard Micro Guard MP Microporous Shoe Cover with Elastic Closure



**Laboratory Instrument India**