



Email : sales@laboratoryinstrumentindia.com

Product Name :
International Splash 1 Chemical Splash Sleeves with Taped
Seams and Elastic Closures

Product Code :
LBNY-00017-11000685



Description :

High moisture vapor transfer rate prevents the buildup of heat and moisture to keep wearers cooler
Possible applications include manufacturing, emergency response, environmental cleanup, and more.
Disposable, chemical splash sleeves are designed with red taped seams to protect your arms from chemical splashes and toxic particles.
Proprietary, hybrid membrane fabric helps keep arms cool as you work
Secure, elastic closures on both ends help keep sleeves in place even with repetitive movements
Longer length protects more of your arm from exposure to workplace hazards
Taped seams offer an extra layer of protection by preventing contaminants from entering or seeping through exposed seams
Elastic closures on both ends of sleeves keep them in place to prevent skin exposure.

Technical Specification :

International Splash 1 Chemical Splash Sleeves with Taped Seams and Elastic Closures



Laboratory Instrument India

