



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
VR Grommet Aprons

**Product Code :**  
LBNY-00017-11000687



**Description :**

Blue aprons feature seamless durability and strength. Ideal for critical and controlled environments. Comfortable, adjustable ties for neck and waist allow for flexible positioning and custom fit. Strong, texturized nylon ties designed for extended wear; interchangeable for better sanitation. Constructed of proprietary VR material from top to bottom, including the grommets, for an environmentally safe, more durable product.

**Technical Specification :**

VR Grommet Aprons



**Laboratory Instrument India**