



Email : sales@laboratoryinstrumentindia.com

Product Name :
Eagle Polyurethane on Nylon Suits: Overalls

Product Code :
LBNY-00017-11000691



Description :

Made of 200 denier nylon
Interior polyurethane coating and exterior durable water repellency coating for 100% waterproof protection
Mildew resistant
Resistant to most hydrocarbon oils, organic acids, salts, alkalies and some organic solvents
Breathable polyurethane on nylon garments for all-day comfort and wear
Tingley's polyurethane on nylon suits are tough and generously cut to fit over the heaviest of work clothes.
Transportation, meat packing, utilities, municipalities and general industry.

Technical Specification :

Eagle Polyurethane on Nylon Suits: Overalls



Laboratory Instrument India