



Email : sales@laboratoryinstrumentindia.com

Product Name :
Spiral Spring

Product Code :
LBNY-0003-1010028



Description :

Spiral Spring

Technical Specification :

Material: Stainless Steel
Usage: Industrial
Surface Treatment: Color Coated
Material Grade: SS304
Technique: Hot Rolled
Overall length: 55mm.
Coiled length: 20mm.
Diameter: 15.5mm
Wire diameter: 0.7mm
Coil opens at 20g with an extension of 40mm for 100g load



Laboratory Instrument India

