



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mens Lightweight Khaki Flame Resistant Cargo Pant

Product Code :
LBNY-00017-11000698



Description :

Mens Lightweight Khaki Flame Resistant Cargo Pant has six total pockets, including a cargo pocket on each leg to ensure you have plenty of room to store all your tools.

Ideal for storage and breathability

Left cargo pocket has two snap closures and right utility pocket is open

Wash Care: Industrial Laundry - Light Soil, Home Wash

Fabric: Flame resistant, 7 oz. (235 g/m²)

Hazard: Flash Fire, Arc Flash

Waistband with hook and eye closure

Blend: 48% Modacrylic/37% Lyocell/15% Aramid

Sits at the waist

Regular fit in thigh and seat with slightly tapered leg

Two front slack-style pockets and two cargo utility pockets

Product Collection: Cool touch 2.

Technical Specification :

Mens Lightweight Khaki Flame Resistant Cargo Pant

Laboratory Instrument India

