



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Light Blue Lab Coat with Exterior Pocket, Grippers

**Product Code :**

LBNY-00017-11000699



**Description :**

The left chest pocket with pencil stall and two lower pockets will keep you stocked with pens and small tools while the one-piece, lined collar and button-front closure will keep you looking professional during rounds and time in the lab.

Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash

Closure: Five grippers

Collar: One-piece, lined, notched lapel

Pocket: Left chest exterior pocket, Bartacked pencil stall, Two lower exterior pockets

Primary Closure: Button

Side vent openings for easy access to underneath attire

Fit Details: Fits true to size. Meant to be worn over clothes.

Non-yellowing buttons for long-lasting whiteness and professional appearance

Requires minimum ironing thanks to a wrinkle-resistant finish.

**Technical Specification :**

Light Blue Lab Coat with Exterior Pocket, Grippers



**Laboratory Instrument India**

