



Email : sales@laboratoryinstrumentindia.com

Product Name :
Superior Glove Contender Aramid Cut-Resistant Sleeves

Product Code :
LBNY-00017-11000702



Description :

Material: Aramid Fiber

Cuff Style: Knit

For Use With (Application): Metal handling, automotive, Metal fabrication, welding

Construction: Aramid Encased, Silica infused Fiber outer layer with Cool max Inner layer

Size: X-Large

Length (English): 22 in.

Includes: Thumb hole.

Technical Specification :

Superior Glove Contender Aramid Cut-Resistant Sleeves



Laboratory Instrument India