



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mens Navy Lapel Counter Coat

**Product Code :**  
LBNY-00017-11000703



**Description :**

Mens poly/cotton lapel counter coat with exterior pockets.

Wash Care : Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash

Fabric: 5 oz. Performance Blend Poplin

Blend: 80% Polyester / 20% Combed Cotton

Finish: Wrinkle Resistant

Closure: Three buttons

Collar: One-piece, lined, notched lapel

Lab coat with plenty of pocket space between a left chest pocket and two lower ones so you can keep tools handy. The men's three button front closure on this coat features non-yellowing buttons.

Non-yellowing buttons for long-lasting whiteness and professional appearance

Requires minimum ironing thanks to a wrinkle-resistant finish

Blazer styling provides a professional look

One chest pocket and two lower pockets provide storage for tools and personal items

Materials: Poplin, Cotton/Polyester.

Primary Closure: Button

Length: 30 in. (based on size medium - back)

Pocket: Left exterior chest pocket, Bartacked pencil stall, Two lower exterior pockets

Industrial laundry friendly

Fit Details: Fits true to size. Meant to be worn over clothes.

Gender: Male

**Technical Specification :**

---

Mens Navy Lapel Counter Coat



**Laboratory Instrument India**