



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
White Nylon Hairnet

**Product Code :**  
LBNY-00017-11000704



**Description :**

Strong, elastic headband  
Prevent stray hairs from contaminating food or work surfaces  
Cool, economical hair coverage option when full containment is not necessary.  
Disposable  
Provides full hair coverage.

**Technical Specification :**

White Nylon Hairnet



**Laboratory Instrument India**