



Email : sales@laboratoryinstrumentindia.com

Product Name :
Long Sleeve Comfort Woven Lightweight Shirt Regular

Product Code :
LBNY-00017-11000707



Description :

Flame resistant long sleeve comfort lightweight regular shirt
Fabric: Bulwark Exclusive 5.5 oz. 50% Aramid 25% Lyocell 25% Modacrylic- Flame-resistant
Topstitched cuffs with two-button adjustable closure
Tailored sleeve placket
Care: Industrial Laundry - Light Soil, Home Wash
Banded, Topstitched collar with red trim at neck sea
Larger patch chest pockets
Hemmed front placket with button closure
Full side seam gusset.

Technical Specification :

Long Sleeve Comfort Woven Lightweight Shirt Regular



Laboratory Instrument India