



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Comfort Lightweight Pant Length: 34 in.

**Product Code :**  
LBNY-00017-11000708



**Description :**

Bulwark flame resistant lightweight pants  
One-piece waistband with concealed button closure  
Two scoop front pockets, concealed pocket on right pocket bag, two patch hip pockets  
Full crotch gusset  
Fabric: Ripstop Flame-resistant, 6.9 oz. (235 g/m<sup>2</sup>) 65% Cotton/ 35% Polyester  
Double utility pocket on right leg  
Double articulated knees one in. silver trim sewn just below reinforced knee  
Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash  
Protection: Arc Rating ATPV 9.0 calories/cm<sup>2</sup>.  
Utility loop extension on front belt loop  
Cargo pocket with snap closures on left leg.

**Technical Specification :**

Comfort Lightweight Pant Length: 34 in.



**Laboratory Instrument India**

