



Email : sales@laboratoryinstrumentindia.com

Product Name :
Comfort Knit Long Sleeve Henley Regular

Product Code :
LBNY-00017-11000709



Description :

Flame resistant long sleeve knit shirt.
Concealed rib knit cuffs
Full underarm gusset
Concealed chest pocket with pencil stall
Concealed pencil stall sleeve pocket
Four-button placket
Rib knit henley collar
Raglan sleeves
Fabric: Exclusively with fabrics by Milliken Flame-resistant, 5.3 oz. (180 g/m²) 69% Cotton / 25% Polyester / 6% Polyoxadiazole
Care: Home Wash
Contrasting thread color seam tag Exterior CAT 2 2112 labels
Protection: Arc Rating ATPV 8.2 calories/cm².

Technical Specification :

Comfort Knit Long Sleeve Henley Regular

Laboratory Instrument India

