



Email : sales@laboratoryinstrumentindia.com

Product Name :
Long Sleeve Comfort Woven Lightweight Shirt Tall

Product Code :
LBNY-00017-11000710



Description :

Fabric: Bulwark Exclusive 5.5 oz. 50% Aramid 25% Lyocell 25% Modacrylic- Flame-resistant
Care: Industrial Laundry - Light Soil, Home Wash
Flame resistant long sleeve comfort lightweight regular shirt
Tailored sleeve placket
Hemmed front placket with button closure
Larger patch chest pockets
Banded, Topstitched collar with red trim at neck seam
Topstitched cuffs with two-button adjustable closures
Full side seam gusset
Seam tag.

Technical Specification :

Long Sleeve Comfort Woven Lightweight Shirt Tall



Laboratory Instrument India

