



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Long Sleeve Comfort Woven Lightweight Shirt Tall

**Product Code :**  
LBNY-00017-11000710



**Description :**

Fabric: Bulwark Exclusive 5.5 oz. 50% Aramid 25% Lyocell 25% Modacrylic- Flame-resistant  
Care: Industrial Laundry - Light Soil, Home Wash  
Flame resistant long sleeve comfort lightweight regular shirt  
Tailored sleeve placket  
Hemmed front placket with button closure  
Larger patch chest pockets  
Banded, Topstitched collar with red trim at neck seam  
Topstitched cuffs with two-button adjustable closures  
Full side seam gusset  
Seam tag.

**Technical Specification :**

Long Sleeve Comfort Woven Lightweight Shirt Tall



**Laboratory Instrument India**

