



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Straight Fit Jean with Stretch Length: 34 in.

**Product Code :**

LBNY-00017-11000711



**Description :**

Flame resistant stretch jeans 34 in. in length

Fabric with stretch

Fabric: Flame-resistant, 12 oz. (405 g/m<sup>2</sup>) 94% Cotton/ 5% Polyester/ 1% Spandex

Care: Industrial Laundry -Light Soil, Home Wash

Felled inseam and seat seam

Soft, pre-washed and sanded fabric

One-piece waistband with metal button closure, metal is not exposed to skin

Five jean-style pockets

Bulwark waistband ticket and size strip

Protection: Arc Rating ATPV 15 calories/cm<sup>2</sup>.

**Technical Specification :**

Straight Fit Jean with Stretch Length: 34 in.



**Laboratory Instrument India**

