



Email : sales@laboratoryinstrumentindia.com

Product Name :
Chill-Its Cooling Arm Sleeves

Product Code :
LBNY-00017-11000715



Description :

Chill-Its Cooling Arm Sleeves provide full arm coverage to absorb and wick away sweat, deliver advanced cooling relief, shield skin from UV rays, and even protect against light abrasions.

Protection from light abrasions

Reflective accents

Anti-odor treatment

Performance knit fabric is soft, stretchable and durable

Wicks moisture when dry to absorb sweat from the body

Instant cooling when wet that can last for hours

UPF 50+ Protection

Machine washable

Chemical-free cooling.

Technical Specification :

Chill-Its Cooling Arm Sleeves



Laboratory Instrument India

