



Email : sales@laboratoryinstrumentindia.com

Product Name :
Chill-Its Phase Change Cooling Vest Elastic Extenders

Product Code :
LBNY-00017-11000716



Description :

Made of 100% elastic with hook and loop extenders
Chill-Its Phase Change Cooling Vest Extenders provide 3.5 in. (9 cm) of extension when used with the Chill-Its Phase Change Vest.
Comfortable stretch
For use with cooling vests.

Technical Specification :

Chill-Its Phase Change Cooling Vest Elastic Extenders



Laboratory Instrument India