



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Chill-Its Phase Change Cooling Vest Elastic Extenders

**Product Code :**  
LBNY-00017-11000716



**Description :**

Made of 100% elastic with hook and loop extenders  
Chill-Its Phase Change Cooling Vest Extenders provide 3.5 in. (9 cm) of extension when used with the Chill-Its Phase Change Vest.  
Comfortable stretch  
For use with cooling vests.

**Technical Specification :**

Chill-Its Phase Change Cooling Vest Elastic Extenders



**Laboratory Instrument India**