



Email : sales@laboratoryinstrumentindia.com

Product Name :
Protective Wear Integrity Cleanroom Coveralls

Product Code :
LBNY-00017-11000727



Description :

Reusable coveralls with ESD qualities provide maximal fluid resistance and particle filtration for a variety of cleanroom applications including pharmaceutical, biotechnology, and medical device manufacturing.
Fluid and bacteria resistant Integrity high performance carbon striped fabric
Full-length, lint-limiting front zipper
Three male snaps on the back of the leg for boot attachment
Mandarin collar with two snaps at top
Laser-cut, pre-serged seams
Silicon-free
ESD carbon stripes provide electrostatic dissipation.

Technical Specification :

Protective Wear Integrity Cleanroom Coveralls



Laboratory Instrument India