



Email : sales@laboratoryinstrumentindia.com

Product Name :
Enterprises Eco Beard/Hair Net Dispenser, Clear Recycled
Plastic with Hinged Lid

Product Code :
LBNY-00017-11000729



Description :

High quality dispenser.
Durable PETG hinges support a lid that keeps contents clean and dust free.
Recycled materials support sustainability initiatives and LEED Credits.
Dispenser uses clear plastic and includes pre-drilled holes for easy mounting.
ZING Beard/Hair Net Dispenser conveniently stores beard or hair nets for quick and easy access.
Beard/Hair Net Dispenser is made from heavy-duty 100% recycled PETG providing superior durability and crack resistance over other materials like acrylic.
Dispenser is ideal for high use applications and can be mounted on a wall or placed on a counter-top.
Measures 9 in.H x 7.5 in.W x 8 in.D. Beard/Hair.

Technical Specification :

Enterprises Eco Beard/Hair Net Dispenser, Clear Recycled Plastic with Hinged Lid



Laboratory Instrument India

