



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Voltmeters

**Product Code :**  
LBNY-0003-1010032



**Description :**

Voltmeters

**Technical Specification :**

Voltmeter are used in the laboratories of schools & Colleges for doing the practicals in Physics lab.  
Case: Abs Plastic Yellow case  
Accuracy: +/- 2% F.S.D  
Desk Stand: Front terminal ABS Plastic Yellow stand  
Double scale, 50mm long to read 0 – 5 V x 0.2 V, in plastic case.



**Laboratory Instrument India**