



Email : sales@laboratoryinstrumentindia.com

Product Name :
Safe-Sense Bouffant Caps

Product Code :
LBNY-00017-11000733



Description :

Bouffants are manufactured from a strong, lightweight, and breathable spun-bounded polypropylene (SPP) fabric to protect against dirt, grime, and certain dry particulates in controlled environments.

21 in. diameter

Choose from pleated or traditional style.

Bouffant caps

Used for cleanroom

Non-sterile

Technical Specification :

Safe-Sense Bouffant Caps



Laboratory Instrument India