



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Safe-Sense Bouffant Caps

**Product Code :**  
LBNY-00017-11000733



**Description :**

Bouffants are manufactured from a strong, lightweight, and breathable spun-bounded polypropylene (SPP) fabric to protect against dirt, grime, and certain dry particulates in controlled environments.

21 in. diameter

Choose from pleated or traditional style.

Bouffant caps

Used for cleanroom

Non-sterile

**Technical Specification :**

Safe-Sense Bouffant Caps



**Laboratory Instrument India**